

CUSTOMER INFORMATION		SOURCE INFORMATION		
COMPANY	DIVISION	CUSTOMER REFERENCE #	NEXT TEST DUE	SEQUENCE #
PURCHASE ORDER	LICENSE	ADDITIONAL SOURCE DESCRIPTION		
CONTACT PERSON		SOURCE MANUFACTURER		
Voice:	Fax:	SOURCE MODEL	SOURCE SERIAL NUMBER	
ADDRESS		SOURCE ASSAY DATE	LEAK TEST INTERVAL	
			days	
		ISOTOPE(S)	ACTIVITY	
			mCi	
		PERSON COLLECTING SAMPLE	SAMPLE COLLECTION DATE	



State of Wisconsin License 071-1234-01

LIMITATION OF LIABILITY

While RAM Services, Inc., maintains comprehensive quality control and data verification procedures, we make no warranty of any kind, expressed or implied, for the consequences of erroneous test results. Neither RAM Services, Inc., nor its employees or agents shall be liable under any claim, charge, demand, whether in contract, tort, or otherwise, for any and all loss, cost, charge, or demand, fee, expense or damage of any nature or kind arising out of, connected with, resulting from, or sustained as a result of any laboratory test requested.

© RAM Services, Inc. 1999

EMERGENCY PROCEDURES

If, during sample collection activities, you believe the integrity of the source has been compromised, take the following actions:

- ☛ MOVE to low-background area.
- ☛ STOP leak-test collection.
- ☛ WARN other personnel in immediate area.
- ☛ SECURE area from unauthorized entry.
- ☛ NOTIFY RSO or RAM Services, Inc.

RAM Services use only.

Analysis Date _____

Sample Activity _____ **µCi**

Analyzed by _____

SEALED RADIOACTIVE SOURCE LEAK-TEST INSTRUCTIONS

Leak Test Kit Model GK-1

- X Obtain Leak-Test Kit and ensure cotton swab, re-sealable plastic bag are enclosed.
- X Obtain a properly functioning, calibrated survey instrument. (If available)
- X Complete or correct the customer information portion of leak-test kit.
- X Put on latex gloves, if available, and, if applicable, dosimetry.
- X Insure source is in the off or shielded position. (If survey meter is available, confirm expected dose rates.)
- X Wipe the location of the source according to RSO/manufacturer instructions. (For Ni-63 gas chromatograph, wipe the air-vent port or location recommended by manufacturer.
- X After completion of wipe, place cotton swab in plastic bag and seal.
- X Survey hands and equipment if meter is available.
- X Remove latex gloves and wash hands.
- X Complete wipe test information portion of this leak-test kit and place in outer shipping package.
- X If survey meter is available, survey shipping package. If exposure rate exceeds 0.5 mR/hr, contact RSO or RAM Services, Inc.

RETURN COMPLETED LEAK TEST KIT TO:

RAM SERVICES INC.

510 COUNTY HIGHWAY V

TWO RIVERS, WISCONSIN 54241

VOICE: +1.920.686.3889

FAX: +1.920.686.3899

LEAKTESTS@RAMSERVICESINC.COM